



Counterlung Size

Two sizes available - Small and Standard

They both have the same lung volume but the Small is designed for narrower shoulders and a shorter back and is designed for use with the Small BMCL harness only.

In order to fit different body sizes the counterlung position can be adjusted vertically by using the alternative mounting holes. It is essential to keep the counterlungs against your back and upper rear shoulder to avoid excessive hydrostatic imbalance.

Feet		Height												
		Under	5'	5' 2.5"	5' 5"	5' 7.5"	5' 10"	5' 12.5"	6' 1"	6' 3.5"	6' 6"	6' 8.5"	Over	
Inches	cm		60	62.5	65	67.5	70	72.5	75	77.5	80	82.5		
			152	159	165	171	178	184	191	197	203	210		
Chest	Under	Small					Standard							
	36	91	Small					Standard						
	38	97	Small					Standard						
	40	102	Small					Standard						
	42	107	Small					Standard						
	44	112	Small					Standard						
	46	117	Small					Standard						
	48	122	Small					Standard						
	50	127	Small					Standard						
	52	132	Small					Standard						
54	137	Small					Standard							
56	142	Small					Standard							
Over		Small					Standard							